



Advanced Test Equipment Rentals
www.atecorp.com 800-404-ATEC (2832)

INSTRUCTIONS

OLYMPUS SCP-10

OES INSTANT MEDICAL CAMERA



OLYMPUS

1 FEATURES AND SPECIFICATIONS

1-1 Features

Mounted on the Olympus fiberscope, the SCP-10 performs instant photography.
The camera features:

- ① On-the-spot recording of your findings.
- ② Compatibility with both OES and non-OES fiberscopes by the use of adapters.
- ③ Standard and magnified images with the OES fiberscope.
- ④ Center-weighted light measurement giving priority to the most important subject area for correct exposure.
- ⑤ LED (light-emitting diode) indicators incorporated on the camera and inside the viewfinder turn on to warn of impossible exposure when the safety cover is ejected and when all the films have been exposed.
- ⑥ An LED indicator on the adapter turns on to verify correct connections between the components – OES fiberscope-adapter-camera – when the electrical connection is completed.

1-2 Specifications

Applicable film	Polaroid High Speed Color Land film, Type 779
Number of exposures per film pack	10
Shutter	Eveready type
Shutter speeds	1, 1/2, 1/4, 1/8, 1/15, 1/30, 1/60, 1/125 sec. and B
Film transport	Motorized, fully automatic transport
Exposure measurement	TTL center-weighted metering
Film counter	Subtracting type; automatic reset
Viewfinder	Single-lens reflex camera with diopter adjustment ring and eyecup
Connection to adapter	Exclusive mount with automatic lock
Connection to fiberscope	OES bayonet mount via the adapter
Shutter release	Push-button type. Cable release is attachable.
Interchangeable adapters	Two adapters (2.6X and 3.6X) for OES fiberscopes, and one for non-OES fiberscopes